

# Fast 2D azimuthal integration using PyFAI



| The European Synchrotron

Jérôme Kieffer, online data analysis @ ESRF

- **X-ray scattering/diffraction experiments using 2D detectors:**
    - X-ray scattering by electrons
    - Benefits of using 2D detectors
  - **Azimuthal integration in python**
    - How azimuthal integration is performed in Python
    - What is PyFAI: installation, testing and basic use
    - Improvements provided by pyFAI
      - **Pixel splitting algorithms**
      - **Parallelization of the algorithms**
  - **Calibration of the experimental setup**
    - Peak picking methods & geometry optimization setup
  - **Extra tools**
- by Giannis Ashiotis
- by Aurore Deschildre


# X-ray scattering and 2D detectors

- Scattering by an electron:**

- $\mathbf{q} = \mathbf{k}_1 - \mathbf{k}_0$

- $A = b \exp(-2i\pi (\mathbf{k}_1 - \mathbf{k}_0) \cdot \mathbf{r} / \lambda)$

where b is the scattering cross-section of an electron (Thomson or elastic scattering)



The cross section of an atom or scattering density is increasing with Z, the number of electrons

- Intensity of scattering: Only the norm of the scattered intensity is accessible:**



- $I(\mathbf{q}) = AA^* = \sum_j \sum_k b_j b_k \exp(4\pi(\mathbf{k}_1 - \mathbf{k}_0) \cdot (\mathbf{r}_j - \mathbf{r}_k) / \lambda)$

- *Sum over all pairs of electrons*

# Small angle X-ray diffraction

## Experimental setup:

- $q = |\mathbf{k}_1 - \mathbf{k}_0| = 4\pi \sin(2\theta/2) / \lambda$
- Measure  $I = f(q)$
- Small deviation:  $2\theta < 5^\circ$
- *Sample detector distance: few meters*
- *Detector orthogonal to incident beam (tilt neglected)*




Data from BM29

Azimuthal integration




PyFAI

Or other  
FIT2D  
SPD  
XRDUA  
Foxtrot ...



# X-ray diffraction and powder diffraction

- **Atom plans in a crystal acts like a grating:**
  - Constructive interference for:
    - $n \lambda = 2d \sin(\theta)$
    - Known as **Bragg's law** (1913)
- **The angle between incident and diffracted beam is discrete, depending on the d-spacing inside the crystal**
  - Angle of diffraction:  $2\theta = 2 \sin^{-1}(n\lambda / 2d)$
- **Powder are composed from randomly oriented cristallites: only those in Bragg condition diffracts**
  - Diffraction along cones centered on the sample
  - Intersection with the detector plan: conic
    - Called **Debye-Scherrer rings**
    - PyFAI is compatible with ellipses, parabola or hyperbola



## Pros:

- Collect more photons thanks to larger solid angle
- Great statistic

## Cons:

- CCD & CMOS detector need dark-current subtraction
- Most detector need flat-field correction to equalize the pixel response
- Optical fiber taper used together with CCD detectors induce distortion
- Each pixel has a different solid angle when viewed from sample
- Polarization of the X-ray beam induces variation of signal
- Intensity can need to be linearized (power  $p \approx 1.0$ )
- Pixel are always too large and limit your resolution.

## PyFAI made the choice of:

- Pixel-wise correction 
$$I_{corr} = \frac{I^p - I_{dark}^p}{flat * \delta\Omega * Polarization * Absorption}$$
- Geometry transformation directly to destination space (r,  $\chi$ ) or (q,  $\chi$ ) or (2 $\theta$ ,  $\chi$ )
  - **This includes the taper distortion: no taper corrected image is produced**

# PyFAI for azimuthal integration



- **Primarily a Python library to perform azimuthal integration**
  - The main object is an AzimuthalIntegrator instance which has:
    - **The geometry of the experimental setup**
    - **High performance integration/averaging routines**
    - **Pixel-wise correction of the image for:**  
dark current, flat field, solid angle, polarization, mask, absorption...
- **Tools to fit the geometry using a reference sample and diffracted rings**
- **Set of detectors (almost 40 different detectors)**
- **Set of reference sample (about 10 different calibrants)**
- **Tools to perform distortion correction for a fibre optic taper**
- **Tools to write easily plugins for other applications**
  - Lima, EDNA, PyMca, ...

- **PyFAI is an open source project ...**
  - Hosted on <http://github.com/kif/pyFAI>
  - Bug and issue tracker on github:
  - Mailing list: [pyFAI@esrf.fr](mailto:pyFAI@esrf.fr)
    - **Open outside ESRF and archived**
- **Dependencies:**
  - Requested
    - **Python 2.6 or 2.7**
    - **Numpy, scipy and matplotlib**
    - **FabIO to be able to read images**
  - Optional
    - **Cython and a C compiler to compile binary modules**
    - **PyOpenCL and a supported GPU for best performances**
    - **Sphinx to build the documentation**

# Tutorial 1: Test pyFAI

- Download PyFAI and test it (under a UNIX computer):

```
$ git clone https://github.com/kif/pyFAI
```

```
$ cd pyFAI
```

```
$ python setup.py build
```

```
$ python setup.py test
```

```
$ python setup.py build_doc
```

```
$ firefox ./build/sphinx/html/index.html
```


```
$ cd benchmark; ./benchmark.py -c -g
```

At ESRF:


```
$ export https_proxy=http://proxy.esrf.fr:3128/
```

```
$ export http_proxy=http://proxy.esrf.fr:3128/
```

← takes a variable amount of time due to downloads



← takes a lot of time but nice graphics



- **Why Python:**
  - Widely used general-purpose, high-level programming language
  - Design philosophy emphasizes code readability: fewer lines of code
  - Clear programs on both a small and large scale.
- **PyFAI code analysis (line of code):**
  - Python 12k lines of code
  - Cython 8k lines of code → generated 400k lines of C
  - OpenCL 2k lines of kernel code → to run on GPU for example
  - Tests: 3k lines of python
- **Most of the tutorial will use Python**
  - For scripting: ipython notebook with the pylab interface:  
*\$ ipython notebook --pylab*

# A bit theory about azimuthal integration

- **Re-grid data on a regular, polar grid**

- Input data:  $x, y$
- Output data:  $(r, \chi)$  or  $(q, \chi)$  or  $(2\theta, \chi)$


$$r = \sqrt{x^2 + y^2}$$

$$2\theta = \tan^{-1}(r/d)$$

$$q = 4\pi \sin(2\theta / 2) / \lambda$$

$$\chi = \tan^{-1}(y/x)$$


**Not linear functions:**  
each needs a specific re-gridding  
or interpolation



# Tutorial: Generate a test image with rings

- **Create a test image with concentric rings**
  - **Work in radius, math are easier but the method is the same as q**
$$r = \text{sqrt}(x^2 + y^2)$$
  - **Image shape: 1000x1000**
  - **center: 500,500,**
  - **rings spaces by 50 pixels**
  - **Gaussian shape for the rings:  $I = \exp(-(r-r_0)^2)$**

# Correction: Generate a test image with rings



IP[y]: Notebook

Untitled0

Last saved: May 07 4:32 PM

File Edit View Insert Cell Kernel Help

Code

```
In [4]: shape = 1000,1000      # shape of image
        xc,yc = 500,500        # position of center
        rings = range(0,500,50) #position of rings in pixel
```


```
In [5]: y,x = numpy.ogrid[:shape[0],:shape[1]]
        x -= xc
        y -= yc
        r = numpy.sqrt(x*x+y*y)
        img = numpy.zeros(shape)
        for radius in rings:
            img += numpy.exp(-(r-radius)**2) #gaussian profile
```

```
In [6]: imshow(img)
```

```
Out[6]: <matplotlib.image.AxesImage at 0x4cb6d50>
```

# Tutorial: Azimuthal integration, first implementation

- **Average the intensity of all pixels with radius in  $[r, r+dr]$** 
  - Get pixel coordinates of all pixel between  $r$  and  $r+dr$  (numpy.where)
  - Take intensities and average them (mean method of numpy array)
  - Loop over  $r$
- **Make a function taking an image and the radius array as input**
  - Should output the  $I = f(\text{radius})$
  - Benchmark it (use %timeit)





# Correction: Azimuthal integration, first implementation

- **One possible solution:**

```
In [41]: def integrate0(img, radius, bins):
        rmin = radius.min()
        rmax = radius.max()
        delta_r = (rmax-rmin)/(2.0*bins)
        integrated = numpy.zeros(bins)
        positions = linspace(rmin+delta_r, rmax-delta_r, bins)
        for i,r0 in enumerate(positions):
            valid = abs(radius-r0) < delta_r
            pix_pos = numpy.where(valid)
            intensities = img[pix_pos]
            integrated[i] = intensities.mean()
        return positions, integrated
```


If you decompose mean = sum / number  
- number can be obtained from histogram  
- sum can be obtained from weighted histogram  
→ numpy has a faster implementation for histograms

```
In [43]: %timeit R,I = integrate0(img,r, 500)
```

1 loops, best of 3: 6.11 s per loop

```
In [42]: clf()
        R,I = integrate0(img,r, 500)
        plot(R,I)
```

```
Out[42]: [<matplotlib.lines.Line2D at 0x840d590>]
```



- **Drawback: Terrible speed !**

- **Look at the documentation of `numpy.histogram`**
- **Perform the histogram of the radius array**
  - What is the meaning of it ?
- **Perform the histogram of the radius array weighted by image intensity**
  - What is the meaning of it ?
  - The mean being the sum divided by the number of element ...  
can you perform the integration ?
- **Re-implement the former function using histograms**
- **What are the performances ?**

# Correction: histogram based integration

```
In [82]: def integratel(img, radius, bins):  
         u_hist, edges = numpy.histogram(radius, bins)  
         w_hist, edges = numpy.histogram(radius, bins, weights=img)  
         positions = (edges[:-1]+edges[1:])/2.0  
         integrated = w_hist/(1.0*u_hist)  
         integrated[u_hist==0] = 0  
         return positions, integrated
```


- **Performances are much better**

```
In [83]: %timeit R,I = integratel(img,r, 500)  
  
1 loops, best of 3: 118 ms per loop
```

```
In [84]: clf();  
         R0,I0 = integrate0_corr(img,r, 500)  
         plot(R0,I0,label="first")  
         R1,I1 = integratel(img,r, 500)  
         plot(R1,I1,label="hist")  
         legend()
```

```
Out[84]: <matplotlib.legend.Legend at 0x98e69d0>
```

- **Note the small differences:**  
**Related to the number of pixels/bin**



## About the precision of this method ...

- **Perform the same analysis when the original image has 20% noise**
- **Perform the analysis data with 100, 200, 500, 1000 and 2000 bins**
  - Overlay all results on a single graph
  - What do you see ?

Two major steps are performed during azimuthal integration:

- Pixel wise corrections
- Regridding on a polar grid

Two main methods:

- *`AzimuthalIntegrator.integrate1d`*
- *`AzimuthalIntegrator.integrate2d`*

Exercise: look at the function signature

# Pre-processing: pixel-wise corrections

- **Masked pixels are not treated ...**
  - Hence mask should be defined on the raw image and not on the distortion corrected image.
    - **Define the mask using the “drawMask\_pymca” script provided**
    - **Perform the “un-correction” of a mask defined for FIT2d:**
      - Use the `pyFAI.distortion.Distortion.uncorrect(img)` method (slow)
- **Intensity correction:**
  - Dark current subtraction:  $I - I_{\text{dark}}$
  - Flat field correction  $I / I_{\text{flat}}$
  - Solid angle correction:  $I / \cos^3(\alpha)$  where  $\alpha$  is the incidence angle
  - Polarization correction:  $I * 2 / ((1 + \cos(2\theta) - f \cos(2\chi) * (1 - \cos(2\theta))))$
  - Linearity correction and absorption correction shall be addressed upstream
- **All calculation are performed in floating point:**
  - single precision with error compensation on GPU (Kahan summation\_
- **All correction can be activated or disable.**
  - Only solid angle correction is activated by default

# Regridding is not interpolation

- **Interpolations are not well adapted to experimental needs:**
  - Most common are nearest neighbor, bi-linear, bi-cubic, ...
  - No conservation of the intensity, nor locally nor globally
- **Integration provides intensity conservation.**
- **Various patterns of pixel splitting implemented in pyFAI:**
  - No pixel splitting (like EMBL Hamburg):
    - **Histogram type algorithms: noisy**
  - Fixed splitting on 2 bins (like Foxtrot, Soleil):
    - **Bins and image size are linked, all pixels have the same size**
  - Bounding box splitting (like FIT2D, ESRF)
    - **Pixels are assumed to be parallel to the output grid**
    - **Over estimation of the pixel size, tend to smear-out features**
  - Full pixel splitting
    - **Pixels are split according to their actual area, assuming they are quadrilateral**

Nota: PyFAI integration algorithms provide intensity conservation only without pixel-wise preprocessing,  
i.e. solid angle correction deactivated

# Tutorial: influence of the pixel splitting scheme

- **Use PyFAI.AzimuthalIntegrator to generate the geometry:**
  - Define a Fairchild detector (4096x4096 detector, 15 $\mu$ m square pixel size)
  - Place the beam-center at (2000, 2000), the detector at 100 mm from sample
  - Generate the radius array:  $r=f(x,y)$
- **Create an image with rings on it**
  - Add noise to the signal
- **Perform integration using no splitting, bounding box and pixel splitting**
- **Look at the result where the statistics are low**
- **Look at the peak shape**



# Correction: influence of the pixel splitting scheme

```
In [6]: import pyFAI, numpy
ai = pyFAI.AzimuthalIntegrator(detector="fairchild") # Fairchild is a 4k x 4k detector with 15 microns pixel size
shape = ai.detector.max_shape # retrieve the size of the image
bins = 2000 # number of bins for integration
ai.setFit2D(100, 2000, 2000) # Place the detector at 100mm with the beam center at 2000, 2000
r = ai.rArray(shape) # contains the radius array
```




```
In [9]: img = numpy.zeros(shape)
sigma = 50e-6 #that >3 pixels
for radius in numpy.arange(0, 0.045, 0.005):
    img += numpy.exp(-(r-radius)**2/(2*sigma**2))
img_noise = img + numpy.random.random(shape) #ratio S/N = 1
imshow(img_noise)
```

```
Out[9]: <matplotlib.image.AxesImage at 0x7f7edd8cc490>
```

```
In [10]: clf();
plot(*ai.integrate1d(img_noise, bins, unit="r_mm",method="numpy"),label="histogram")
plot(*ai.integrate1d(img_noise, bins, unit="r_mm",method="splitBBox"),label="splitBBox")
plot(*ai.integrate1d(img_noise, bins, unit="r_mm",method="splitPixel"),label="splitPixel")
legend()
```


```
Out[10]: <matplotlib.legend.Legend at 0x7f7eddb6a990>
```

```
In [11]: #difference between splitBBox and splitPixel:
R,I_bbox = ai.integrate1d(img_noise, bins, unit="r_mm", method="splitBBox")
R,I_pixel = ai.integrate1d(img_noise, bins, unit="r_mm", method="splitPixel")
plot(R,I_pixel-I_bbox)
```



# Statistical consideration on pixel splitting ...

- **Influence of pixel splitting on errors:**
  - 1000 frames taken without sample (water only) and average to get
    - $\chi^2$  is too low when using pixel splitting (regardless the scheme)
      - **Correction specific for each geometry (under study)**



Question: Do you know the Point Spread function of your detector ?

# Direct integration(histogram) vs look-up table (SmDv)

- **Histogram or direct integration:**
  - Efficient on a single core (serial execution)
  - Hard to parallelize due to write conflicts
    - **Either have large locked region (atomic operation)**
    - **Either transform the problem is a parallel problem**
  - Well suited to changing geometries
- **“Scatter to Gather” transformation**
  - Store the contribution of each pixel to each bin in a matrix
  - Processing similar to a Sparse Matrix Dense Vector (SmDv) multiplication
- **Look-up table integration or backwards integration:**
  - Easy to parallelize, highly efficient on GPU
  - Need to calculate the LUT in a “serial” way
    - **Not suited if geometry changes often !**

- **Create a detector object for a Titan from Oxford diffraction**
  - See the factory or the class ...
- **Create an AzimuthalIntegrator object with the detector and the geometry**
  - Play with set getFit2D/setFit2D/getPyFAI/setPyFAI
- **Create a calibrant object for Lanthanide hexaboride**
  - Use the fake\_calibration\_image to generate an image
- **Try out integrate1d and integrate2d**
  - Plenty of options, all described in the documentation
  - Benchmark it
- **Issue: how do I get the geometry in real life ???**
  - Use the calibration: next chapter

# Programs to perform integration: pyFAI-integrate

- **pyFAI-integrate: the graphical interface for integration**
- **pyFAI-waxs: command line interface for integration**
- **pyFAI-saxs: command line interface for integration**
- **diff\_tomo: diffraction mapping & tomography tool**

# Graphical user interface to integrate a set of images

1) Click here  
for file dialog  
to load PONI file

Poni File:  ... save to File

Detector:  Wavelength (m):

Pixel1 (m):  Pixel2 (m):

Spline file:  ...

Distance (m):  Rotation 1 (rad):

Poni 1 (m):  Rotation 2 (rad):

Poni 2 (m):  Rotation 3 (rad):

☐ Mask File  ...

☐ Dark Current  ...

☐ Flat Field  ...

☐ Dummy value  delta dummy

☒ Polarization factor  ☒ Solid Angle corrections

Radial units: ☒ 2 $\theta$  ( $^{\circ}$ ) ☐ 2 $\theta$  (rad) ☐ q (1/nm) ☐ q (1/Å) ☐ r (mm)

Number of radial points:  ☐ Std-err (Poisson law)

☐ Number of azimuthal points (2D)  ☐  $\chi$  discontinuity at 0

☐ Radial range

☐ Azimuthal range

☒ Use OpenCL Platform:  Device:

0% Help Reset Save Cancel OK

Parameters  
populated from  
poni-file

Check to  
activate  
correction

File dialogue  
To select images

Mandatory

Chose your  
output space

GPU selection

- **Usage:** `pyFAI-waxs [options] -p ponifile file1.edf file2.edf ...`
- Options:
  - `--version` show program's version number and exit
  - `-h, --help` show this help message and exit
  - `-p PONIFILE` PyFAI parameter file (.poni)
  - `-n NPT` Number of points in radial dimension
  - `-w WAVELENGTH, --wavelength=WAVELENGTH`  
wavelength of the X-Ray beam in Angstrom
  - `-e ENERGY, --energy=ENERGY`  
energy of the X-Ray beam in keV ( $hc=12.398419292\text{keV.A}$ )
  - `-u DUMMY, --dummy=DUMMY`  
dummy value for dead pixels
  - `-U DELTA_DUMMY, --delta_dummy=DELTA_DUMMY`  
delta dummy value
  - `-m MASK, --mask=MASK` name of the file containing the mask image
  - `-d DARK, --dark=DARK` name of the file containing the dark current
  - `-f FLAT, --flat=FLAT` name of the file containing the flat field
  - `-P POLARIZATION_FACTOR, --polarization=POLARIZATION_FACTOR`  
Polarization factor, from -1 (vertical) to +1  
(horizontal), default is None  
for no correction, synchrotrons are around 0.95
  - `--error-model=ERROR_MODEL`  
Error model to use. Currently on 'poisson' is  
implemented
  - `--unit=UNIT` unit for the radial dimension: can be  $q_{\text{nm}}^{-1}$ ,  $q_{\text{A}}^{-1}$ ,  
 $2\theta_{\text{deg}}$ ,  $2\theta_{\text{rad}}$  or  $r_{\text{mm}}$
  - `--ext=EXT` extension of the regrouped filename (.xy)

- **usage: pyFAI-saxs [options] -n 1000 -p ponifile file1.edf file2.edf ...**

Azimuthal integration for SAXS users.

- positional arguments:

FILE                List of files to calibrate

- optional arguments:

-h, --help           show this help message and exit

-v, --version

-p PONIFILE          PyFAI parameter file (.poni)

-n NPT               Number of points in radial dimension

-w WAVELENGTH, --wavelength WAVELENGTH  
                     wavelength of the X-Ray beam in Angstrom

-e ENERGY, --energy ENERGY  
                     energy of the X-Ray beam in keV (hc=12.398419292keV.A)

-u DUMMY, --dummy DUMMY  
                     dummy value for dead pixels

-U DELTA\_DUMMY, --delta\_dummy DELTA\_DUMMY  
                     delta dummy value

-m MASK, --mask MASK name of the file containing the mask image

-d DARK, --dark DARK name of the file containing the dark current

-f FLAT, --flat FLAT name of the file containing the flat field

-P POLARIZATION\_FACTOR, --polarization POLARIZATION\_FACTOR  
                     Polarization factor, from -1 (vertical) to +1  
                     (horizontal), default is None for no correction,  
                     synchrotrons are around 0.95

--error-model ERROR\_MODEL  
                     Error model to use. Currently on 'poisson' is  
                     implemented

--unit UNIT          unit for the radial dimension: can be q\_nm^-1, q\_A^-1,  
                     2th\_deg, 2th\_rad or r\_mm

--ext EXT            extension of the regrouped filename (.dat)

pyFAI-saxs is the SAXS script of pyFAI that allows data reduction (azimuthal integration) for Small Angle Scattering with output axis in q space.



# Integration tool: Diffraction mapping/tomography

## Tool: `diff_tomo --help`

### Description of the options:


```
diff_tomo      #Name of the program
-o dt.h5       # output file
-t 153         # number of translation
-r 181         # number of rotations
-d dark.edf    # image with dark-current
-p Si.poni     # File with the calibration
-g            # Perform the calculation on the GPU
data/*.edf     # all 2D diffraction images
```

→ see tutorial

Performances measured on the processing: 100 ms/frame

40Mpix/s → overhead of the IO !!!

Reading/writing 217GB over an hour → 75MB/s



# Geometry determination

Come back to real world:  
Given a reference sample 2D diffraction pattern,  
determine the geometry of the experiment


- **The determination of the geometry is also known as calibration**
- **PyFAI assumes this setup does not change during the experiment**
- **It is divided into 4 major steps:**
  - Pre-processing of images: averaging, dark/flat correction, desaturation
  - Identification of peaks and groups of peaks belonging to same ring
  - Least-squares refinement of the geometry parameters on peak position
  - Validation by an human being of the geometry

# Pre-processing tool: PyFAI-average

- **Average out many files using various functions:**
  - Mean: by default
  - Median: remove zinger
  - Max/min → see example
  - Any other function available from ndarray stacks
- **Average all but outliers**
  - cutoff based on standard deviation

Do this for your:

- Data image
- Dark-current image
- Flat field image



- **Flat-field images integrate the signal response of each pixel**

→ **It is not the image of the beam without sample !!!**

- They are acquired using an isotropic source (fluorescence)
- Multiple detector/source configuration should be averaged
- The flat image is rather stable in time but varies with energy !
- The flat field image should ideally have a average signal of 1


- **PyFAI can estimate the contribution from the taper:**

- Taken into account in the solid-angle correction
- Needs to be deactivated if you provide a flat-field image

```
AzimuthalIntegrator.set_correct_solid_angle_for_spline(False)
```

# Pre-processing tool: masking tool

- **Most detectors have some invalid areas:**
  - Inter-module space (Pilatus detectors)
  - CCD pixels not connected to the taper
    - **pyFAI provides a set of 39 detectors (68 labels) with their masks**
- **You will need to mask the beamstop + other shadows**
- **drawMask\_pymca**
  - Provided by PyFAI
  - Requires PyMca4 installed



# Note about masks ... when using a distorted detector

- **PyFAI can use FIT2D masks but ...**
  - FIT2D masks are drawn on distorted corrected images
  - PyFAI applies the mask on the raw image
- **If you have a distorted detector and a mask drawn with FIT2D:**
  - Either draw a new mask, using `drawMask_pymca`
  - Or distort your mask to apply it on the raw image:

```
import pyFAI, fabio, pyFAI.distortion
fit2d_msk = fabio.open("fit2d.msk").data
detector = pyFAI.detectors.FReLoN("eo2k.spline")
distortion = pyFAI.distortion.Distortion(detector)
pyfai_mask = distortion.uncorrect(fit2d_msk)[0]>0
edf = fabio.edfimage.edfimage(data=pyfai_mask.astype('int8'))
edf.save('pyfai_mask.edf')
```


- **Retrieve the default mask for a detector:**

```
import pyFAI, fabio
det = pyFAI.detectors.detector_factory("mar345").astype('int8')
fabio.edfimage.edfimage(data=det.mask).write('mar345_mask.edf')
```

# Geometry optimization: least squares refinement

- **Geometry:**

- 3 distances, in meter: dist, poni1, poni2
- 3 angles in radians: rot1, rot2, rot3 (later not changed usually)
- Optionally the wavelength can be optimized (correlated to dist at low  $2\theta$ )



- **Optimization:**

- Optimization of square of the distance in  $2\theta$  to the calibrant in radians
  - **Can weight optimization using intensities**
- Minimize the cost function using Sequential Least Squares Programming
  - **`scipy.optimize.fmin_slsqp`**



# Peak picking: extract peak position from image

- **Steepest ascent**
  - Starting from a seed, search for the nearest (local) maximum
  - Sub-pixel precision obtained from second order expansion:  $\delta_{pos} = -H^{-1} \times \nabla$ 
    - **Methods to calculate the hessian and gradient, noise issue**
- **Massif algorithm**
  - Groups of peaks are extracted from a difference with a blurred image
  - Unique parameter: width of the gaussian for blurring
  - Uses binning and un-binning to keep processing time reasonable
- **Blob-detection: systematic peak search**
  - Difference of subsequent blurs. Search for all
  - Octave changing with binning → see Aurore's presentation
- **Monte-Carlo sampling**
  - Peak selection when region is known
  - Based on random search + steepest ascent

# Calibrants: provide opening of cones

- **PyFAI ships 10 reference sample (decreasing  $2\theta$  of first ring):**
  - Au: Gold
  - CeO<sub>2</sub>: Ceria
  - Si: Silicon
  - $\alpha$ -Al<sub>2</sub>O<sub>3</sub>: Corundum
  - Cr<sub>2</sub>O<sub>3</sub>: Chromium oxide
  - CrO<sub>x</sub>: Undefined chromium oxide formerly used on MX beamlines
  - LaB<sub>6</sub>: Lanthanide hexaboride
  - PBBA: Para Bromo Benzoic Acid
  - C<sub>14</sub>H<sub>30</sub>O: tetradecanol
  - AgBh: Silver Behenate
- **But you can provide your d-spacing file if you prefer:**
  - Ascii text files with d-spacing written in Angstrom (like FIT2D)
  - Use the American Mineralogist database:
    - <http://rruff.geo.arizona.edu/AMS/amcsd.php>

# Strength of the geometry

- **Can accommodate highly tilted detector**
  - See tutorial on LaB6
- **Can import/export geometry of FIT2D**
  - getFit2D/setFit2D method from geometry/azimuthalIntegrator
- **Stable to binning:**
  - LaB6 at 0.1m from detector (60 $\mu$ m, 2048\*2048), orthogonal geometry
  - Calibration performed after binning: precision  $\sim 1/10$  of a pixel

Binning	Pixel Size ( $\mu$ m)	Poni1 (m)	Poni2 (m)	Delta1 ( $\mu$ m)	Delta2 ( $\mu$ m)
1	60	6.00E-002	6.00E-002	4.820	1.630
2	120	6.00E-002	6.00E-002	-8.530	-4.020
4	240	6.00E-002	6.00E-002	-2.700	-4.210
8	480	6.00E-002	6.01E-002	45.410	91.440

- **Few examples from tests:**
  - HalfCCD FReLoN image of LaB6 from ID11 taken at  $\lambda=0.27\text{\AA}$
  - FReLoN image of Si from ID22 taken at  $\lambda=0.42\text{\AA}$
- **Extra tutorial: See in folder with examples**
  - Calibration of images taken with the detector on a 2 $\theta$ -arm
    - **LaB6 images taken with the Titan detector at  $\lambda=1\text{\AA}$**
  - Saxs data with a pixel detector (dealing with gaps in detector)
    - **AgBh images taken with a Pilatus1M at  $\lambda=1\text{\AA}$**

# Additional tools:

Extra scripts which can make your life easier

# List of scripts in pyFAI:

- **Pre-processing tools:**
  - drawMask\_pymca: tool for drawing a mask
  - pyFAI-average: tool for averaging/median/... filtering images
- **Calibration tools:**
  - pyFAI-calib: Initial calibration tool
  - pyFAI-recalib: Calibration refinement tool (Obsolete)
  - MX-calibrate: Calibrate automatically a set of images
  - check\_calib: checks the calibration of an image at the sub-pixel level
- **Azimuthal integration tools:**
  - pyFAI-integrate: graphical interface for integration
  - pyFAI-waxs: command line interface for integration
  - pyFAI-saxs: command line interface for integration
  - diff\_tomo: diffraction mapping&tomography tool

# Calibration tool: MX-calibrate

- **Tool to calibrate a set of diffraction images taken at various distances (especially MX beamlines)**

- Expect distance to be in the filename for linear regression
- Can read the wavelength from file header (or from command line)


- **Example of use:**

*MX-calibrate -D Pilatus6M -c CeO2 ceria\_\*.cbf*

- **Returns: Linear regression of parameters vs detector distance**

```
dist = 0.00100143591695 * dist_mm + -0.00112366879462    R= 0.999999965812    stderr= 1.17109288577e-07
poni1 = -3.41861282886e-06 * dist_mm + 0.215648741605    R= -0.996810176908    stderr= 1.2240640473e-07
poni2 = 1.64908246321e-06 * dist_mm + 0.212896001186    R= 0.935403434948    stderr= 2.78771824164e-07
rot1 = 1.17701119402e-06 * dist_mm + 0.00689408012249    R= 0.676514144642    stderr= 5.72993860345e-07
rot2 = -1.56845501269e-06 * dist_mm + 0.00502610867189    R= -0.895473163114    stderr= 3.48664084642e-07
rot3 = 3.11858790559e-12 * dist_mm + -1.13933614082e-09    R= 0.60959568452    stderr= 1.81362213889e-12
direct = 1.00147388349 * dist_mm + -1.12406896453    R= 0.99999996682    stderr= 0.000115373737414
tilt = 9.14175326104e-06 * dist_mm + 0.487856144005    R= 0.118853078555    stderr= 3.41542499598e-05
trp = 0.0128891131751 * dist_mm + 143.975004565    R= 0.936685043656    stderr= 0.00215489991677
centerX = -0.0350625303913 * dist_mm + 1238.48253482    R= -0.99989808711    stderr= 0.000223882943102
centerY = 0.00427419793753 * dist_mm + 1254.36314238    R= 0.975374520063    stderr= 0.000432230388786
```

- **Tool to validate the quality of a calibration:**
  - Perform the azimuthal integration
  - Project back in the averaged pattern on the original image
  - Shows the difference map, highlighting:
    - **Distortion spline file errors**
    - **Defects in the dark measurment**
    - **Inhomogeneities in the detector response (flat)**
  - Calculate the offset between the two images using phase correlation
    - **Precision better then 0.1pixel (here  $\delta x = +0.07$  &  $\delta y = -0.001$ )**



→ proper mask is mandatory



- **PyFAI can deal with distortion correction for detectors:**
  - Spline-file generated from FIT2D
  - Displacement maps for Pilatus detectors (in percent of pixel)
- **They are obtained from a regular grid**
  - There is no tool (yet) to generate displacement maps in pyFAI
  - Can be added if needed → 2 weeks of work
- **PyFAI can re-sample data on a regular grid if needed:**
  - Similar to azimuthal integration (full pixel splitting)

```
import pyFAI, fabio, pyFAI.distortion
det = pyFAI.detectors.FReLoN('splinefile.spline')
dis = pyFAI.distortion.Distortion(det)
cor_img = dis.correct(raw_img)
```

- **EDNA:**
  - Generic plugin to perform azimuthal integration
  - Optimized plugin for BM29 BioSaxs (10 Mpixel/s in production)
- **LlmA and ProcessLib :**
  - Azimuthal integration performed online (→ BM01)
  - Distortion correction
  - Demo available ...
- **PyMca**
  - Under construction
- **Plugins based on a worker (pyFAI.worker.Worker instance):**
  - Azimuthal integrator + information about input and output shapes
  - Contains dark/flat/...
  - Communicates using JSON strings
  - Should be initialized for best performances
  - Output using a writer: pyFAI.io.Writer instance offering EDF, HDF5, ...

## Appendix

- **Coding standards**
  - Comments: every class/method has a documentation
  - Keep the code and the API simple for developers & users
- **Code development**
  - Most functions are internally tested to prevent regression
    - **However we are lacking functional tests on scripts**
  - Continuous integration:
    - **run test daily to prevent regression**
      - **The *master* branch is always useable**
- **Question about versionning:**
  - Would you like version number like 2014.06 ?
  - Version 1.0 is likely to never arrive

# Local installation and bootstrapping

- **Test are always run from a local install of the library obtained from:**

```
$ python setup.py build
```

- **They are located in build/lib\*/pyFAI**

- **One can use this build to test programs thanks to F.Picca:**

```
$ ./bootstrap ipython
```

```
In [1]: import pyFAI
```

```
In [2]: pyFAI
```

```
Out[2]: <module 'pyFAI' from '/users/kieffer/workspace-  
ssd/pyFAI/build/lib.linux-x86_64-2.6/pyFAI/__init__.py'>
```

- **1200 commits since it started in July 2011**
  - 11 contributors from 5 synchrotrons world-wide + academic labs
    - **ESRF, Soleil, Petra3, APS, Sesame + CEA, ...**
- **Languages used:**
  - Python 23k lines
    - **including 8k of Cython → generates 400.000 lines of C**
    - **Including 3k for the tests**
  - OpenCL: 2k lines for parallel kernels
- **Mantra: K.I.S.S**
  - Simple is easier to understand
  - Simple is easier to maintains
  - Try to be self contained: easier to deploy !

# Thanks to all contributors:

- **ESRF**
  - Jonathan Wright
  - Manuel Sanchez del Rio
  - V. Armando Solé
  - Aurore Deschildre
  - Amund Hov
  - Peter Boesecke
- **LinkSCEEM:**
  - Dimitris Karkoulis
  - Zubair Nawaz
  - Giannis Ashiotis
- **Synchrotron Soleil:**
  - Frédéric Emmanuel Picca
- **Advanced Photon Source:**
  - Clemens Prescher
- **Petra3:**
  - Gunthard Benecke
  - Gero Flucke



| The European Synchrotron